

Dental Needle

Product Information Card



Description:

Sterile dental needle for single use. The dental needle is used as an injector needle of cartridge syringes for dental for dental anaesthesia and other purposes. Does not contain components made of natural rubber latex.

Components:

Description	Raw Materials
Cap	Polypropylene (+ 5% Polyethylene)
Hard Case	Polypropylene
Hub	Polyethylene + master batch
Cannula	Stainless steel
Lubricant	Silicone
Adhesive	Epoxy resin



Product Range:

Product Code Imperial Thread <i>(with the hub)</i>	Product Code Metric Thread <i>(blue hub)</i>	Colour	Cannula Ø		Cannula Length	
					Effective	At cartridge
	DN*2508F	Red	25G	0.50 mm	8 mm	10.7 mm
	DN*2535F	Red	25G	0.50 mm	35 mm	10.7 mm
	DN*2716F	Yellow	27G	0.40 mm	16 mm	10.7 mm
DN*2721		Yellow	27G	0.40 mm	21 mm	10.7 mm
	DN*2721F	Yellow	27G	0.40 mm	21 mm	10.7 mm
	DN*2721B	Yellow	27G	0.40 mm	21 mm	13.2 mm
DN*2722		Yellow	27G	0.40 mm	22 mm	10.7 mm
DN*2730		Yellow	27G	0.40 mm	30 mm	10.7 mm
DN*2735		Yellow	27G	0.40 mm	35 mm	10.7 mm
	DN*2735B	Yellow	27G	0.40 mm	35 mm	13.2 mm
DN*2741		Yellow	27G	0.40 mm	41 mm	10.7 mm
	DN*3013F	Green	30G	0.30 mm	13 mm	10.7 mm
	DN*3013B	Green	30G	0.30 mm	13 mm	13.2 mm
DN*3016		Green	30G	0.30 mm	16 mm	10.7 mm
	DN*3016F	Green	30G	0.30 mm	16 mm	10.7 mm
	DN*3016B	Green	30G	0.30 mm	16 mm	13.2 mm
DN*3019		Green	30G	0.30 mm	19 mm	10.7 mm
DN*3021		Green	30G	0.30 mm	21 mm	10.7 mm
	DN*3021F	Green	30G	0.30 mm	21 mm	10.7 mm
	DN*3021B	Green	30G	0.30 mm	21 mm	13.2 mm
DN*3022		Green	30G	0.30 mm	22 mm	10.7 mm

Dental Needle

Packaging: **Single unit:** Cap/case and sealing label
Unit box: 100 pieces
Shipping carton: 2000 pieces
Dimensions =
334 x 264 x 213 mm for items of 13.2 mm cannula at cartridge
284 x 264 x 213 mm for items of 10.7 mm cannula at cartridge

Sterilization: With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Storage: Shelf life of 5 years
* The batch number and expiry date are printed on the sealing label
* Do not store at extreme temperature and humidity

Precautions for Use: Check integrity of the unit package before use
For single use only
Do not bent the needle prior to use, as breakage and possible patient-injury may result
Dispose of safely after single use to avoid risk of infection

Status of Regulation: Class IIa, sterile product following MDD 93/42/EEC

Notified body: TÜV Rheinland LGA Products, GmbH 0197
Manufacturer: TERUMO CORPORATION JAPAN
2-44-1 Hatagaya, Shibuya-ku
Tokyo, 151-0072 Japan
EC representative: TERUMO EUROPE N.V.
Interleuvenlaan 40
B - 3001 Leuven
Belgium



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Terumo Medical Products EMEA
+32 16 38 12 11

EMEA Sales Offices

Terumo Europe NV African Business Division
+32 16 38 13 08

Terumo Europe NV Benelux Sales Division Belgium
+32 16 339 25 80
The Netherlands
0800 0220396

Terumo Europe NV Emerging Markets Division
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 80 230

Terumo Deutschland GmbH Switzerland
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 130 96 13 00

Terumo Italia S.R.L.
+39 0651 96 14 20

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

Terumo Sweden AB Denmark
+45 7020 93 80

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480 440

Terumo BCT Tibbi Cihazlar Dagim ve Hizmetleri A.S.
+90 216 590 0074

