

**LARIDENT S.R.L.** UNIPERSONALE  
**VIA LAMANEIGRA ,12 R**  
**16030 TRIBOGNA (GE)**  
**TEL.+390185938068 FAX. +390185938595**

# TECNICAL SHEET

## *OF PRODUCT*

|                           |                                                                                                                                                                                        |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>NAME PRODUCT</b>       | PROTHESIS BOX WITH INTERNAL FILTER AND LID                                                                                                                                             |
| <b>CODE</b>               | G1                                                                                                                                                                                     |
| <b>MATERIAL</b>           | POLYPROPYLENE                                                                                                                                                                          |
| <b>THERMAL PROPERTY</b>   | AUTOCLAVABLE UP TO 140°C                                                                                                                                                               |
| <b>PACKING</b>            | IT IS PACKED INDIVIDUALLY IN PLASTIC BAG                                                                                                                                               |
| <b>CHARACTERISTIC</b>     | IT IS AVAILABLE IN THE FOLLOWING COLORS:<br>WHITE, AZURE, FUXIA, YELLOW, GREEN, BLUE .<br>MEASURE: DIAM 9 cm - H 8,2 cm<br>IT RESISTS TO ALL DISINFECTION MATERIALS CURRENTLY IN USE . |
| <b>SAFETY DATA SHEET</b>  | N/A                                                                                                                                                                                    |
| <b>WARNING</b>            | NONE                                                                                                                                                                                   |
| <b>ANALYSIS OF RISK R</b> | N/A                                                                                                                                                                                    |
| <b>CLASSIF. CEE 93/42</b> | NO MEDICAL DEVICE                                                                                                                                                                      |

N/A NO APPLICABLE