

	TECHNICAL DATA SHEET PRIVATE LABEL			
BRAND	check off	Reference code	Ref.CADENCE	Description
HYGITECH		HY-80038		Hygitest Test Prion
DENTO-VIRACTIS				
Name of the product (manufacturer)	Class 6 Emulator - 134C - 18 min			
Manufacturer Reference	HY-80038			
Directive Class (under Directive 93/42/EEC)	Not Applicable		Rule of classification	
Location of manufacture	Istanbul - Turkey			
Products description	Class 6 Emulator - 134C - 18 min, box with 250 Emulators to test Steam sterilisation conditions inside packs.			
Packaging description	Carton box containing the Hygitest Test Prion Emulator strips			
Technical specifications	Refer to technical file			
Composition / ingredients	Paper and Ink			
Instructions for use	Not Applicable			
Security data sheet	Not Applicable			
Technical File	Technical file HY-80038			
Shelf life of the product	Guarantee : 5 years after production date			
Precautions for handling or storage	See Carton / Packaging			
Traceability	Through Lotnumber at the Carton			
Conditioning				
units per CARTON /sale unit	550			
number of cartons per layer of the pallet	4			
number of cartons per pallet	16			Pallet height max : 2,40m
units per pallet	8800			
Supplier	SP Medikal San Ltd Sti			
Contact	Ms. Seda Kucukyilmaz			
Address	Topcular Mah, Incirlik Sokak, Ertaş İş Merkezi 5/ 62-64 Kat 2 Rami/Eyup Istanbul Turkey			
Telephone number	+90 212 613 80 54			
Fax	+90 212 613 80 55			
E-mail Address	seda@spmedikal.com			
Delivery lead time after receipt of subsequent orders	2 weeks			
Delivery lead time (after order reception)	10 weeks			
Transport conditions	Ex-works - Istanbul - Turkey			
Transport costs	Not applicable			
Packaging costs	Included in Product price			
Minimal quantities - Printing (units)	?			
Minimal quantities - Production (units)	?			
Minimal quantities - Order (units)	?			
Payment	100% at despatch date.			

On behalf of HYGITECH
M. Michael Stemmer
Date :

On behalf of
Mr Peter M. den Uil
Date : 31/12/2014


 key solutions for today's world
 SP Medikal San. Ltd. Şti.
 Topcular Mah. Incirlik Sok. Ertaş İş Merkezi No:5/62-64 Kat:2 Rami/Eyüp
 İstanbul-TÜRKİYE Tel: +90-212-613 80 54 Fax: +90-212-613 80 55
 Bayrampaşa V.D.: 781 046 3595 info@spmedikal.com www.spmedikal.com