



<b>Designation</b>	<b>Hydrophilic gauze swabs sterile Medibase</b>
<b>Reference(s)</b>	<b>3135500</b>
<b>Product picture</b>	
<b>Product description</b>	<ul style="list-style-type: none"> <li>• Dimensions folded: 5*5 cm</li> <li>• Number of threads: 17</li> <li>• Number of plies: 12</li> <li>• Sterilisation: Ethylene Oxide</li> </ul>
<b>Characteristics and composition</b>	<ul style="list-style-type: none"> <li>• Dimensions folded: 5*5 cm</li> <li>• Number of threads: 17</li> <li>• Number of plies: 12</li> <li>• Sterilisation: Ethylene Oxide</li> </ul>
<b>Packaging</b>	Carton of 100 swabs (50 sterile pouches of 2 swabs)
<b>Shelf life</b>	5 years
<b>Indications of use</b>	
<b>Legal information</b>	<p>Medical Device of class I sterile under European Directive 93/42/EEC including Directive 2007/47/CE            CE 0123 Marking- Certifying body: TUV PRODUCT SERVICE            Manufacturer: FOMED Industries Inc</p>