



Designation	Sterile non-woven swabs Medibase			
Reference(s)	Size	5x5 cm	5x5 cm	7,5x7,5 cm
	Box of	150 swabs	250 swabs	80 swabs
	Internal conditioning	75 pouches of 2 swabs	50 pouches of 5 swabs	40 pouches of 2 swabs
	Reference	3135502	3135501	3135503
Product picture				
Product description	Non-woven swabs Medibase packed under sterile pouches			
Characteristics and composition	<p>Swabs 60% viscose / 40% polyester</p> <ul style="list-style-type: none"> • Number of plies : 4 • Sterilisation : ethylene oxide 			
Packaging	Available in two sizes (5x5 cm and 7,5x7,5 cm) and three conditioning			
Shelf Life	5 years			
Indications of use				
Legal information	<p>Medical Office of class I unofr European Directive 93/42/EEC including Directive 2007/47/CE CE 0123 Marking. Certifying body : TUV PRODUCT SERVICE Manufacturer: FOMED Industries Inc</p>			