


Designation	Micro-applicators			
References	SIZE	COLOR	PACKAGING	REFERENCE
	REGULAR	BLUE	Box of 4 tubes each containing 100 micro-applicators	3101006
		BLUE/PINK		3101007
		BLUE/ORANGE		3101008
	FINE	BLACK	Tube of 100 micro-applicators	3101009
SUPER FINE	WHITE	3101010		
Product picture				
Product description	Disposable micro-applicators with foldable head.			
Characteristics and composition	<p>Polypropylene handle, with breakable tip: allows to differentiate the micro-Applicators between them!</p> <p>Nylon bristles.</p> <p>Overall length: 102 mm</p> <p>Tube packaging: hygienic dispensing</p>			
Packaging	<p>REGULAR : box of 4 tubes each containing 100 micro-applicators</p> <p>FINE / SUPER FINE : tubes of 100 micro-applicators</p>			
Shelf life	N/A			
Indication of use	<p>Regular fibers, non-absorbent and lint-free. The solutions are held in suspension between the fibers to prevent dripping and allow perfect application. Easy access to all areas of application thanks to the foldable head.</p> <p>Disposable: discard after use.</p> <p>For use with silanes, adhesives, etching gel, hemostatic solutions, cavity bottoms and plaque developers.</p> <p>Store away from light and dust in a cool, dry place.</p>			
Legal informations	Class I medical device under European Directive 93/42 / EEC including Directive 2007/47 / EC CE marking.			

Version 1.0

Updated the 2016/07/26 by SMA