

LARIDENT S.R.L.UNIPERSONALE
VIA LAMANEIGRA, 12 R
161030 TRIBOGNA (GE) ITALIA
TEL.0185938068 FAX.0185938595 info@larident.it

TECNICAL SHEET

OF PRODUCT

NAME PRODUCT	COTTON ROLL DISPENSER
CODE	A4
MATERIAL	ABS PLASTIC , RUBBER GASKKETS, TRANSARENT POLISTIRENE
THERMAL PROPERTY	NO AVAILABLE
INSTRUCTION	STERILIZE BEFORE USING . THE DISPENSER WAS STUDIED TO CONTAIN 60 COTTON ROLLS N. 2 IT CAN BE USED WITH COTTON ROLL N. 1
CHEMICAL CHARACTERISTIC	IT RESISTS TO ALL STERILIZATION LIQUID
SAFETY DATA SHEET	N/A
WARNING	NONE
ANALYSIS OF RISK R	N/A
CLASSIF. CEE 93/42	NO MEDICAL DEVICE

N/A NO APPLICABLE