



Product	Paper/plastic bibs																
References	Color	Blue	Green	White	Orange	Pink	Dark blue	Lime	Lavender								
	Reference	116728	116747	116748	3014611	3014612	3107041	3107043	3107042								
Picture																	
Description	Aprons with very absorbent tissue paper bonded to a layer of plastic, for excellent atiers protection.																
Characteristics & Technical data	<table border="1"> <tr> <td>Dimensions</td> <td>Width: 54 cm Height: 60 cm</td> </tr> <tr> <td>Diameter</td> <td>10,5 cm (± 0,5 cm)</td> </tr> <tr> <td>Weight</td> <td>850 g (± 10%)</td> </tr> <tr> <td>Traceability</td> <td>By lot number printed on the roll</td> </tr> </table>									Dimensions	Width: 54 cm Height: 60 cm	Diameter	10,5 cm (± 0,5 cm)	Weight	850 g (± 10%)	Traceability	By lot number printed on the roll
Dimensions	Width: 54 cm Height: 60 cm																
Diameter	10,5 cm (± 0,5 cm)																
Weight	850 g (± 10%)																
Traceability	By lot number printed on the roll																
Storage	Roll of 80 pre cut aprons with tie. Packaged in clear plastic with non-adhesive label. Keep in a cool, dry place.																
Garantee	Not applicable																
Shelf life																	
Legal Information	Not applicable																