

	Technical file	Version: <b>2020</b>	Creation date: 08.02.2020
	<b>Patient Apron No. 3</b>		Last Editing (Nr.): 7
page 1 of 1			

<b>product name</b>	THIENEL Patient Apron, plastic with tissue
<b>reference(s)</b>	30060 - white, 30160 - light-blue, 30260 - light-green, 30360 - dark-orange, 30460 - orange, 30560 - pink, 30660 - yellow, 30760 - dark-blue, 30860 - sky-blue, 30960 - lilac, 31060 – burgundy, 31260 – grey, 31360 - lime
<b>description</b>	disposable apron to protect clothing; tie fastening
<b>composition</b>	polyethylene (9 µm), tissue (20 g/m <sup>2</sup> )
<b>size / length x width (tolerance)</b>	60 x 54 cm
<b>classification</b>	none
<b>instructions for use</b>	1. Tear off apron from roll at perforation, 2. Release ties along perforation, 3. Tear off neckline, 4. Unfold apron and put it on, 5. Fasten ties with a knot behind the neck
<b>unit</b>	roll of 80 pieces
<b>weight unit (tolerance)</b>	850 g (+/- 10 %)
<b>packaging unit</b>	wrapped in transparent plastic foil; leaflet
<b>volume sales unit</b>	0,0792 m <sup>3</sup> *
<b>sales unit</b>	box of 20 rolls
<b>size sales unit / height x length x width</b>	55 x 30 x 48 cm
<b>gross weight sales unit</b>	18,5 kg
<b>storage conditions</b>	keep dry and cool
<b>minimum durability</b>	3 years

\* average volume of product (unit) = volume sales unit : number of units