


Frame + Protection Shields

Technical data sheet

1. Administrative information about the company		Date of issue : 01/10/2015
1.1	Name : Hygitech SAS	
1.2	Complete address : 25 rue Bleue 75009 PARIS – FRANCE	Phone: +33 1 48 01 32 89 Fax : +33 1 48 01 32 80 e-mail : info@hygitech.fr Web site : www.hygitech.fr
1.3	Contact information Materiovigilance : Claudia Régnier	Phone : +33 1 45 91 31 82 Fax : +33 1 48 65 67 23 e-mail : qualite@hygitech.fr
2. Information about the equipment		
2.1	Trade name: Frame + Protection Shields	
2.2	Class of MD : N/A European Directive: 93/42/EEC including Directive 2007/47/CE According Annex n° N/A Notified body: N/A Date of first launch on the UE market: 2015_02 Manufacturer : Hygitech	
2.3	<p>Product description :</p> <p>Complete article: Autoclavable frame, single use shields</p> <ul style="list-style-type: none"> Anti-fog treatment Vertical inclination Frame available in 3 colours: white, blue, pink <p>Box with 1 frame Box of 12 shields</p> <ul style="list-style-type: none"> Comfortable frame, it adapts to all face types Handy: Shields easy to fit Perfect protection: 30cm shield protects the whole face The shield assures clear vision, Abrasion and scratch-resistant <p>Light: 22 grams</p>	
		

2.4	Packaging : <table border="1" style="width: 100%; border-collapse: collapse; margin-bottom: 10px;"> <thead> <tr style="background-color: #f4a460;"> <th style="width: 25%;">Reference</th> <th style="width: 25%;">OU Order unit</th> <th style="width: 25%;">CARTON</th> <th style="width: 25%;">MOQ Minimum order quantity</th> </tr> </thead> <tbody> <tr> <td>HY- 07001</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> </tr> <tr> <td>HY- 07002</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> </tr> <tr> <td>HY- 07003</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> <td>Box with 1 frame</td> </tr> <tr> <td>HY- 07009</td> <td>Box of 12 shields</td> <td>Box of 12 shields</td> <td>Box of 12 shields</td> </tr> </tbody> </table> Labelling : Unknown	Reference	OU Order unit	CARTON	MOQ Minimum order quantity	HY- 07001	Box with 1 frame	Box with 1 frame	Box with 1 frame	HY- 07002	Box with 1 frame	Box with 1 frame	Box with 1 frame	HY- 07003	Box with 1 frame	Box with 1 frame	Box with 1 frame	HY- 07009	Box of 12 shields	Box of 12 shields	Box of 12 shields							
Reference	OU Order unit	CARTON	MOQ Minimum order quantity																									
HY- 07001	Box with 1 frame	Box with 1 frame	Box with 1 frame																									
HY- 07002	Box with 1 frame	Box with 1 frame	Box with 1 frame																									
HY- 07003	Box with 1 frame	Box with 1 frame	Box with 1 frame																									
HY- 07009	Box of 12 shields	Box of 12 shields	Box of 12 shields																									
2.8	Composition of the product, dimensions and accessories : <table border="1" style="width: 100%; border-collapse: collapse; margin-bottom: 10px;"> <thead> <tr style="background-color: #f4a460;"> <th style="width: 50%;">Elements</th> <th style="width: 25%;">Materials</th> <th style="width: 25%;">Dimensions</th> </tr> </thead> <tbody> <tr> <td>HY-07001 – White frame</td> <td>Polypropylene</td> <td>22 grams</td> </tr> <tr> <td>HY-07002 - Blue frame</td> <td>Polypropylene</td> <td>22 grams</td> </tr> <tr> <td>HY-07003 - Pink frame</td> <td>Polypropylene</td> <td>22 grams</td> </tr> <tr> <td>HY-07009 - 12 refill shields</td> <td>Polypropylene</td> <td>30cm shield</td> </tr> </tbody> </table> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr style="background-color: #f4a460;"> <th colspan="3" style="text-align: center;">Presence of :</th> </tr> </thead> <tbody> <tr> <td style="width: 60%;">Latex</td> <td style="width: 20%; text-align: center;"><input type="checkbox"/> YES</td> <td style="width: 20%; text-align: center;"><input checked="" type="checkbox"/> NO</td> </tr> <tr> <td>Phthalates (DEHP)</td> <td style="text-align: center;"><input type="checkbox"/> YES</td> <td style="text-align: center;"><input checked="" type="checkbox"/> NO</td> </tr> <tr> <td>Products of animal or biological origin</td> <td style="text-align: center;"><input type="checkbox"/> YES</td> <td style="text-align: center;"><input checked="" type="checkbox"/> NO</td> </tr> </tbody> </table>	Elements	Materials	Dimensions	HY-07001 – White frame	Polypropylene	22 grams	HY-07002 - Blue frame	Polypropylene	22 grams	HY-07003 - Pink frame	Polypropylene	22 grams	HY-07009 - 12 refill shields	Polypropylene	30cm shield	Presence of :			Latex	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO	Phthalates (DEHP)	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO	Products of animal or biological origin	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO
Elements	Materials	Dimensions																										
HY-07001 – White frame	Polypropylene	22 grams																										
HY-07002 - Blue frame	Polypropylene	22 grams																										
HY-07003 - Pink frame	Polypropylene	22 grams																										
HY-07009 - 12 refill shields	Polypropylene	30cm shield																										
Presence of :																												
Latex	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																										
Phthalates (DEHP)	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																										
Products of animal or biological origin	<input type="checkbox"/> YES	<input checked="" type="checkbox"/> NO																										
2.9	Field - Indications: Unknown																											
3. Method of sterilizing																												
	Sterile : NO Sterilization method: N/A																											
4. Storage conditions																												
	Normal conditions of storage: N/A Special precautions: N/A Shelf life: N/A Presence of temperature indicators: N/A																											
5. Safety in use																												
5.1	Technical security: N/A																											
5.2	Biosecurity: N/A																											
6. Advices of use																												
	Refer to the manual for more information on the instructions for use.																											
7. Additional product information																												
	Single use: YES (shields only)																											
8. Related documents																												
	Unknown																											