

## 1175 SAFETOUCH® ADVANCED™ SLIM NITRILE POWDER-FREE GLOVE

Class I, Medical Device and Category III, Personal Protective Equipment

### FEATURES & BENEFITS

- Great elasticity
- Good dexterity
- Its smooth appearance gives it correct touch and grip of objects
- Made from a thinner nitrile material, Slim Nitrile technology provides great protection and flexible “latex-like” comfort
- Watertight and airtight



### Recommendation:

- Glove for professional use
- In the healthcare and industrial environment, this glove is used to make barrier against microorganisms (bacteria and mold) between the wearer of the glove and the examined people or objects / materials handled.
- Protective glove designed to protect the user from the chemicals tested.
- Powder-free, it reduces cross contamination and limits the risk of irritation.
- For food handling / preparations (HACCP)

### PRODUCT DESCRIPTION

Product name :	<b>SafeTouch® Advanced™ Slim</b>
Product type :	Single use, non-sterile
Material :	Nitrile
Colour :	Blue
Shape :	Ambidextrous
Cuff :	Rolled
Surface :	Textured fingertips
Powdered :	No
Palm thickness :	0,08 mm (± 0,02 mm)
Finger thickness :	0,10 mm (± 0,03 mm)
Unit weight (M):	4 g (± 10%)
Country of origin :	Asia
AQL :	< 1,5
Strength at break:	
✓ Before aging	≥6N
✓ After aging	≥6N

Article	Taille	Longueur	
1175-A	T 6	XS	
1175-B	T 7	S	≥ 240
1175-C	T 8	M	≥ 240
1175-D	T 9	L	≥ 250
1175-E	T 10	XL	≥ 260

### INDUSTRIES

Health / Medical / Dentist / Tattoo / Food / Non residential catering / Community / Administration / Sanitary / Maintenance / Cleaning / Lab / Pharmaceutical / Research / Industry



# 1175 SAFETOUCH® ADVANCED™ SLIM NITRILE POWDER-FREE GLOVE

Class I, Medical Device and Category III, Personal Protective Equipment

## CERTIFICATION

In accordance with requirements of medical devices Dir.93/42/EEC modified by Dir.2007/47/EC

Comply with the EN 455-1/2/3/4 Medical gloves for single use

Manufactured under ISO 13485 & ISO 9001 certified quality system

In accordance with requirements of personal protective equipment Reg. (UE) 2016/425

Comply with the EN 420, EN 374-1/2/4/5(\*) & EN 16523-1 for protective gloves against dangerous chemicals and microorganisms.

In accordance with requirements of Reg.(EC) 1935/2004 and decree of 9th November, 1994. This glove is suitable for food contact(\*\*).

\* Not controlled against viruses.

\*\* Restriction: Do not put in contact with very fatty and acidic foods (see Contact Food Annex - Rev06).

## STORAGE CONDITIONS

Normal conditions of conservation & storage: not to be exposed to moisture and sun, must be stored at a temperature between 5°C and 40°C

Shelf life of the product: 5 years

## LOGISTIC INFORMATION








**Specifications  
outer case**



**Specifications  
inner case**

**Colour code  
inner and outer**

Item	Size mm	Gross weight kg	QTY/ Pallet	Size mm	QTY	Colour
1175-A		4				
1175-B		4,2				
1175-C	325 x 260 x 245	4,5	60 cartons (6 layers of 10 cartons)	235 x 125 x 63	10 boxes of 100 pieces	
1175-D		4,9				
1175-E		5,2				

### EAN Codes

Item	Outer case	Inner case
1175-A	3 700 381 514 062	3 700 381 514 116
1175-B	3 700 381 514 079	3 700 381 514 123
1175-C	3 700 381 514 086	3 700 381 514 130
1175-D	3 700 381 514 093	3 700 381 514 147
1175-E	3 700 381 514 109	3 700 381 514 154