



DATA SHEET 30.215/L

BLACK NITRILE POWDER-FREE GLOVE 240 MM

FIELD OF APPLICATION: industry / chemical / petrochemical / food / away from home food service lab / pharma / research / tattoo / safety field / beauty / jeweller's craft

TECHNICAL SPECIFICATIONS:



Name : HOPEN SAFE BLACK NITRILE GLOVES

Composition: nitrile, powder-free, non-sterile

AQL: 1,5

Colour : black

Texture: finger textured

Form: ambidextrous

Edge: rolled

Tested for food contact: TUV 719195016-MEC11/01-WBH

Origin: Malaysia

Product dimensions:

SIZE		Width mm	Length mm
T 6/7	S	85 ± 5	240 ± 5
T 7/8	M	95 ± 5	240 ± 5
T 8/9	L	105 ± 5	240 ± 5
T 9/10	XL	115 ± 5	240 ± 5



Tensile strength:

Before aging = 6.4 N

After accelerated aging = 6.2 N

PRODUCT SELLING POINTS :

Nitrile material:

- It possess high abrasion, excellent tear and puncture resistance properties, as well as great resistance to broad range of chemicals
- Provide excellent feel, fit, comfort and flexibility to users, resulting in reduced hand fatigue
- Non-allergenic properties of this quality Nitrile glove offer a great alternative to latex

Powder-free glove:

- Ideal for users who are allergic to corn starch powder
- Prevents cross-contamination
- Recommended for protected environment

Ambidextrous glove:

- Can be used on right or left hand

Beaded Cuff:

- The rolled edge of the cuff improves tear glove resistance
- Prevents the cuff from rolling back in use and helps prevent droplets travelling beyond the glove



BLACK NITRILE POWDER-FREE GLOVE 240 MM

REGULATORY:

Class I, CE.

EN 420 : Protective gloves meet general requirements and testing methods

EN 455 - 1/2/3 : Single- use Medical Examination gloves

EN 374 - 2 : Protective gloves against chemicals and micro-organisms

PACKAGING:



Packaging: Quantity: 100 pcs / dispenser box
Type of packaging: carton
Number of dispenser box / carton: 10
Total quantity per carton: 1000 pcs

Carton:

Size of carton (mm): 340 x 270 x 240

Gross weight: 5,3 kg

Net weight: 4,30 kg

Dispenser box:

Size of dispenser box (mm): 230 x 130 x 65

Gross weight: 488 gr

Net weight: 430 gr

Weight of glove pair: 10 gr

Palletization: 54 cartons of 10 dispenser box i.e 54 000 gloves
Pallet weight: 288,4 kg

STORAGE CONDITIONS:

Store in a cool, dry place. Avoid excessive heat (40°C). Protect from ultraviolet sources (sunlight, fluorescent lighting, etc.)

Usage: disposable, single- use

Shelf-life: 5 years

SIMILAR REFERENCES

Reference	Description	QTY
30.202	Blue nitrile powder-free glove - 240 mm	10 x 100
30.201	Blue nitrile slightly powdered glove - 240 mm	10 x 100
30.205	Almond nitrile powder-free glove - 240 mm	10 x 100



Bld de la chanterie - BP 10059 – 49181 Saint Barthélemy d'Anjou Cedex
T : 33 (0) 241 963 434 – F : + 33 (0) 241 963 453
sales@kolmi-hopen.com ; www.h-open.com

