

### Conception et Fabrication de Produits Médicaux et Paramédicaux

## Safety Data Sheet

(REACH regulation (EC) n° 1907/2006-n° 453/2012)  
Revision N°2 Dated 15-02-2017

#### 1. Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

Code:

**GUTTA PERCHA POINTS**

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.3. Details of the supplier of the safety data sheet

Registered company name :

**CFPM**

Address :

2 bis Chemin du loup  
932970 Tremblay en France. France.  
France

Telephone :

+33 (0)1 45 91 31 82

Fax :+(0)1 48 65 67 23

e-mail :

[info@cfpm-medical.fr](mailto:info@cfpm-medical.fr)

##### 1.4. Emergency telephone number:

+ (0) 1 45 42 59 59

INRS / ORFILA <http://www.centres-antipoison.net>

#### 2. Hazards identification.

2.1 Hazardous Ingredient: Zinc Oxide

2.2 PEL/TLV: 5mg/rn' (Zinc fume)

2.3 Dot Hazard Classification: None

2.1 Whmis Classification: None

#### 3. Composition/information on ingredients.

3.1 To main: In usual use, do not present an inhalation, ingestion or contact health hazard.

3.2 To environment: N/A

#### 4. First aid measures.

4.1 Inhalation: N/A

4.2 Skin contact: N/A

4.3 Oral tissue contact: N/A

4.4 Eye contact: consult an ophthalmologist in case of injury

4.5 Ingestion: not hazardous if one or a few are ingested, otherwise consult a physician

#### 5. Firefighting measures.

5.1 Flash Point: None

5.2 Flammable Limits: N/A

5.3 Autoignition Temperature: N/A

5.4 Extinguishing media: Water. Use any means suitable for extinguishing surrounding fire.

5.5 Special Fire Fighting Procedures: Positive pressure, self-contained breathing apparatus.

5.6 Unusual Fire and Explosion Hazards: None

#### 6. Accidental release measures

6.1 Personal precaution: None

6.2 Environmental precaution: None

6.3 Method for cleaning up: Brush

### Conception et Fabrication de Produits Médicaux et Paramédicaux

#### 7. Handling and storage.

**7.1 Handling and Storage:** Keep in closed container stored in a cool, dry and ventilated area. Precautions:  
Avoid exposure to natural light and heat over 38t.

**7.2 Steps to Be Taken if Material is Released or Spilled:** Wear gloves. Pick up spilled materials and place spilled materials  
in container for disposal.

**7.3 Waste Disposal Method:** Follow ail government regulations.

**7.4 Shelf Life:** 4 years from manufacturing date

**7.5 Other Precautions:** Wash hands after use. Store away from food.

#### 8. Exposure controls/personal protection.

**8.1 Respiratory Protection:** Net necessary under normal conditions.

**8.2 Ventilation:** No special ventilation required.

**8.3 Protective Gloves:** Impervious gloves are recommended,

**8.4 Eye Protection:** Chemical safety goggles are recommended.

**8.5 Other Protective Clothing or Equipment:** Common uniform and equipment.

**8.6 Work/Hygienic Practices:** Use normal care when working with chemicals.

#### 9. Physical and chemical properties.

**9.1 Boiling Point:** N/A

**9.2 Vapor Pressure:** N/A

**9.3 Evaporation Rate:** N/A

**9.4 Odor:** None

**9.5 Specific Gravity:** 2.7 g/cm'

**9.6 Vapor Density:** N/A

**9.7 Melting Point:** 100r

**9.8 Solubility in Water:** None

**9.9 Appearance:** Pink Uniformly Conical Shape

#### 10. Stability and reactivity.

**10.1 Stability:** Stable

**10.2 Conditions to Avoid:** Natural light, excessive heat and humidity (Stability)

**10.3 Incompatibility:** Strong oxidizing agents

**10.4 Hazardous Decomposition or By products**

: Carbon Monoxide, Carbon Dioxide, Hydrocarbons and Zinc

: Fumes when heated to decomposition.

**10.5 Hazardous Polymerization:** Will not occur

**10.6 Conditions to Avoid (Polymerization):** None

#### 11. Toxicological information.

**11.1 Not reported**

#### 12. Ecological information.

**12.1 N/A**

#### 13. Disposal considerations.

**13.1 waste receiver**

**14. Transport information.**

**14.1 regulated for shipping:** None

**14.2 identification number:** N/A

**15. Regulatory information.**

**15.1 follow the regulations of the FDA (Food and Drug Administration)**

**16. Other information**

**16.1 For dental professional use only**