




# Technical Data Sheet

Designation	H Files						
References	<b>Length 21 mm</b>						
	N°	08	10	15	20	25	30
	Réf.	281886	281887	281888	281889	281890	281891
	N°	35	40		Assorted- 15-40	Assorted- 45-80	
	Réf.	281892	281893		281894	281895	
	<b>Length 25 mm</b>						
	N°	08	10	15	20	25	30
	Réf.	281896	281897	281898	281899	281900	281901
	N°	35	40		Assorted- 15-40	Assorted- 45-80	
	Réf.	281902	281903		281904	281905	
Visual							
Description	<p>Hedstroem Files, ideal for good root canal shaping</p> <ul style="list-style-type: none"> <li>-Their spiral shape allows very efficient evacuation of debris.</li> <li>-They largely meet the criteria of torsional resistance.</li> </ul> <p>Plastic box of 6 instruments</p>						
Technical Data	<p>Conicity: 2%</p> <p>Non-working blunt tip: this safety tip allows: to ensure smooth guidance during catheterisation, to bypass obstacles without "catching", to avoid incidents and accidents during shaping (no more stoppers, false directions, false canals), to respect: canal curvatures and apical region</p> <p>Materials used:</p> <ul style="list-style-type: none"> <li>- Handle: sterilizable plastic</li> <li>- Stop: silicone</li> <li>- Blade: stainless steel</li> </ul> <p>The instruments are sterilizable at 135°C for 18 minutes.</p>						
Shelf life	<p>Maximum recommended number of uses, under normal conditions:</p> <ul style="list-style-type: none"> <li>- single use for numbers 8, 10 and 15.</li> <li>- 5 for the other numbers</li> </ul>						
Legal information	<p>Class I medical device according to the European Directive 93/42/EEC including Directive 2007/47/EC            Manufacturer: CFPM</p>						

Version 1.0

Mise à jour 22/06/2022 par Nelly S.

**CADENCE – 2 bis chemin du Loup 93290 Tremblay-en-France – FRANCE – Tel. +33 (0)1 49 63 56 20**