



Product Designation	DISPOSABLE SYRINGES FOR SINGLE-USE																		
Reference	<table border="1"> <tr> <td>Content</td> <td>2,5 ML</td> <td>5 ML</td> <td>10 ML</td> <td>20 ML</td> </tr> <tr> <td>Reference</td> <td>3125103</td> <td>3125105</td> <td>3125102</td> <td>3125104</td> </tr> </table>	Content	2,5 ML	5 ML	10 ML	20 ML	Reference	3125103	3125105	3125102	3125104								
Content	2,5 ML	5 ML	10 ML	20 ML															
Reference	3125103	3125105	3125102	3125104															
Product picture																			
Product description	<ul style="list-style-type: none"> ▪ Disposable syringes for single-use for endodontic canal irrigation 																		
Characteristics and composition	<ul style="list-style-type: none"> ▪ Luer-lock thread for a safe fixation of the needle ▪ Latex free piston ▪ Sterile individual package 																		
Packaging	<table border="1"> <thead> <tr> <th>Content</th> <th>Graduation</th> <th>Number of syringes per box</th> </tr> </thead> <tbody> <tr> <td>2.5 ML</td> <td>0.1 ML</td> <td>100</td> </tr> <tr> <td>5 ML</td> <td>0.2 ML</td> <td>100</td> </tr> <tr> <td>10 ML</td> <td>0.5 ML</td> <td>100</td> </tr> <tr> <td>20 ML</td> <td>1 ML</td> <td>80</td> </tr> </tbody> </table>				Content	Graduation	Number of syringes per box	2.5 ML	0.1 ML	100	5 ML	0.2 ML	100	10 ML	0.5 ML	100	20 ML	1 ML	80
Content	Graduation	Number of syringes per box																	
2.5 ML	0.1 ML	100																	
5 ML	0.2 ML	100																	
10 ML	0.5 ML	100																	
20 ML	1 ML	80																	
Shelf life	5 years from the manufactured date																		
Indications of use																			
Legal information	Medical Device of class IIA under European Directive 93/42/EEC including Directive 2007/47/CE CE Marking- 0197-TUV RHEINLAND Manufacturer: CDM CENTER OF EXCELLENCE LTD																		