



# Fiche Technique

Designation	Limes K						
References	<b>Longueur 21 mm</b>						
	N°	06	08	10	15	20	25
	Réf.	0281940	0281941	0281942	0281943	0281944	0281945
	N°	30	35	40	Assortis- 15-40		
	Réf.	0281946	0281947	0281948	0281949		
	<b>Longueur 25 mm</b>						
	N°	06	08	10	15	20	25
	Réf.	0281950	0281951	0281952	0281953	0281954	0281955
	N°	30	35	40	Assortis- 15-40	Assortis 45-80	
	Réf.	0281956	0281957	0281958	0281959	0281960	
Visual							
Description	<p>K-Files: Root canal preparation (marking and reaming). Plastic box of 6 instruments</p>						
Technical Data	<p>Conicity: 2% Non-working blunt tip: this safety tip allows: to ensure smooth guidance during catheterisation, to bypass obstacles without "catching", to avoid incidents and accidents during shaping (no more stoppers, false directions, false canals), to respect: canal curvatures and apical region</p> <p>Materials used: - Handle: sterilizable plastic - Stop: silicone - Blade: stainless steel</p> <p>The instruments are sterilizable at 135°C for 18 minutes.</p>						
Shelf life	<p>Maximum recommended number of uses, under normal conditions: - 2 for the numbers 06 to 20 - 5 for the other numbers</p>						
Legal information	<p>Class I medical device according to the European Directive 93/42/EEC including Directive 2007/47/EC Manufacturer: CFPM</p>						

Version 1.0

Mise à jour 22/06/2022 par Nelly S.