

COMPOSITE & CERAMIC POLISHERS

KENDA NOBILIS 

1-step diamond polishing system


One polisher only for cost effective finishing and polishing of composite restorations and ceramics. The flexibility of the diamond polisher allows a gentle to teeth clean finishing and polishing of composite restorations and ceramics. Perfect smooth surfaces will be achieved.


Recommendation:





KENDA UNICUS 1-Step Diamond Polisher for ultimate high shine polishing (page K-7).



Set REF 0300



Order REF	0301	0302	0303	0304	0305	0306	0307	0308	0309	0310
Size ØxL mm	5,4x13	5x8	10x2	5x8	6,5x10	3x7	3x6	5x10	10x1	9x8,5
ISO No.	808 204...	243 512 055	243 512 050	304 512 100	030 512 050	030 512 060	243 512 030	235 512 035	243 512 050	303 512 100 030 512 090
Grit	fine									
Color	violet									
	7500 – 10000 RPM									

-  Box of: 6, 12, 25, 100
-  Blister of: 6x1, 6x2, 72x1, 72x2
-  Pressure: 0,3 – 2 N
-  Autoclavable: 135°C, 275°F
- Shank: RA, stainless steel
- Medical device: class IIa