




KENDA UNICUS

1-STEP DIAMOND POLISHING SYSTEM

- Polish porcelain and composites with the same polisher
- For an esthetic high luster result
- Ideal for anterior restorations
- No paste for high shine polishing needed
- Recommendation: Use fine grit finishing instruments or KENDA Nobilis before application



SET REF 0350.006A

10



| | | | | | | | | | |
|---------------------------|-------------|-------------|---------------|-------------|---------------------|-------------|-----------------|-------------|-------------|
| REF reorder | 0351 | 0352 | 0353 | 0355 | 0356 | 0358 | 0359 | 0360 | 0361 |
| ISO 1: 838 204... | 243 503 055 | 254 503 050 | 304 503 100 | 030 503 060 | 243 503 030 | 243 503 050 | 303 503 100 | 030 503 090 | 030 503 060 |
| $\varnothing \times L$ mm | 5,4 x 13 | 5 x 8 | 10 x 2 | 6,5 x 10 | 3 x 7 | 5 x 10 | 10 x 1 | 9 x 8,5 | 6 x 7,5 |
| μ grit | extra-fine | | color | | orange | rec. | 5000 - 7500 RPM | max. | 10.000 RPM |
| Refill Box of: | N Pressure: | | Autoclavable: | | Shank: | | Medical device: | | |
| 6, 12, 25, 100 | 0,3 - 2N | | 135°C, 275°F | | RA, stainless steel | | Class IIa | | |

