

KENDA MICROFILL

2-STEP POLISHING SYSTEM

- Optimized for microfilled composites
- Well proven basic polishing system in 2 steps
- Yellow polishers are also ideal to remove cement remnants of orthodontic bands (brackets)
- Available in a great variety of shapes



SET REF 750.012A



15

REF	reorder	4001	4003	4005	4006	4008	4009	4010	4011	5001	5003	5005	5006	5008	5009	5010	5011
ISO	638 204	243 523 055	304 523 100	030 523 060	243 523 030	243 523 050	303 523 100	030 523 090	030 523 060	243 493 055	304 493 100	030 493 060	243 493 030	243 493 050	303 493 100	030 493 090	030 493 060
	∅ × L mm	5,4x13	10x2	6,5x10	3x7	5x10	10x1	9x8,5	6x7,5	5,4x13	10x2	6,5x10	3x7	5x10	10x1	9x8,5	6x7,5
	μ grit	medium		yellow		7500 - 10 000 RPM		10 000 RPM		ultra fine		blue		5000 - 7 500 RPM		10 000 RPM	
	Refill Box of:	N Pressure:		Autoclavable:		Shank:		Medical device:									
	6, 12, 25, 100	0,3 - 2N		135°C, 275°F		RA, stainless steel		Class IIa									

