

LARIDENT SRL

VIA LAMANEIGRA 12 R 16030 TRIBOGNA (GE) ITALY
TEL.00390185938068 FAX. 00390185938595 info@larident.it

TECHNICAL SHEET

OF PRODUCT

PRODUCT NAME	BURS HOLDERS
CODE	LS12 / LS18 / LS24 / LS36 / LS48 / LS72
MATERIAL	BASE: POLYPROPYLENE LOADED MARBLE POWDER COVER: POLYSTYRENE
THERMAL PROPERTIES	THE BASE IS AUTOCLAVABLE UP TO 130 ° C THE COVER IS NOT AUTOCLAVABLE
MECHANICAL, PHYSICAL PROPERTIES OR CHARACTERISTICS	ALL MODELS ARE AVAILABLE WITH BASE FOR HOLES: FG, RA, MIXED (M), ALSO AVAILABLE WITH THE COVER UP (HP). THE BASE RESISTS TO CHEMICAL STERILIZANTS
MEDICAL DEVICES	IT IS NOT A MEDICAL DEVICE BECAUSE THE CONTAINER DOES NOT INFLUENCE THE THERAPEUTIC ABILITY.
PACKING	IT IS PACKED INTO A CARTON BOX
CAUTION	THE COVER IS NOT AUTOCLAVABLE
RISKS ANALYSIS REF.D.L.46 24/2/97	NO
CLASSIFICATION 93/42 EEC	NO MEDICAL DEVICE

TRIBOGNA, 14/02/2017