



Product	CLINIBITE																																		
Reference	FAST: 3000468 EXTRA FAST: 0116244																																		
Picture																																			
Description	<p>Clinibite is an addition curing vinylpolysiloxane especially made for bite registration/</p> <ul style="list-style-type: none"> • Thixotropic, with fast setting time, this high precision material has a special final hardness (95 SHORE A) and is dimensionally stable. Its opacity, its colour and hardness enable the operator to check and use it very easily. 																																		
Characteristics & Technical data	<p><u>Properties:</u></p> <ul style="list-style-type: none"> • High thixotropy: easy control of the application on occlusal surfaces. • Precision: 2 microns: faithfully reproduces the smallest details of occlusion. • High dimensional stability: 0.05% after 24 hours. • Avoids model repositioning errors and overocclusions. • Cutting, easy milling. • Comfortable for the patient. • No distortion <table border="1" data-bbox="349 1626 1485 1939"> <thead> <tr> <th data-bbox="352 1630 472 1749">Technical data Données techniques Datos técnicos Dados técnicos Dati tecnici Technische Daten Technische gegevens</th> <th data-bbox="472 1630 579 1749">Shore A hardness Dureté Shore A Dureza Shore A Durezza Shore A Shore A Härte Hardheid Shore A</th> <th data-bbox="579 1630 735 1749">Details reproduction Reproduction des détails Reproducción de los detalles Reprodução de pormenor Riproduzione dei dettagli Detailwiedergabe Reproductie van details</th> <th data-bbox="735 1630 890 1776">Linear dimensional change Variation dimensionnelle Cambio dimensional Variação dimensional Variazione dimensionale Dimensionsänderung Dimensionale variatie (>24h)</th> <th data-bbox="890 1630 959 1749">Colour Couleur Color Cor Color Farbe Kleur</th> <th data-bbox="959 1630 1094 1749">Mixing ratio Rapport de mélange Relación de mezcla Relação de mistura Rapporto di miscelazione Anmischverhältnis Mengrapport</th> <th data-bbox="1094 1630 1214 1749">Working time * Temps de travail * Tiempo de trabajo * Tempo de trabalho* Tempo de lavorazione* Verarbeitungszeit* Werktijd*</th> <th data-bbox="1214 1630 1342 1749">Setting time in mouth Temps de prise en bouge Tiempo de toma en boca Tempo de presa na boca Tempo di permanenza in bocca Verweildauer im mund Tijd mondopname</th> <th data-bbox="1342 1630 1481 1749">Total setting time* Temps de prise totale* Tiempo de fraguado total* Tempo de presa total* Tempo di presa totale* Gesamtverarbeitungszeit* Totale opname tijd*</th> </tr> </thead> <tbody> <tr> <td data-bbox="352 1776 472 1861">Clinibite Fast</td> <td data-bbox="472 1776 579 1861">95</td> <td data-bbox="579 1776 735 1861">20 µm</td> <td data-bbox="735 1776 890 1861">-0,05%</td> <td data-bbox="890 1776 959 1861">• Green • Vert • Verde • Grün • Groen</td> <td data-bbox="959 1776 1094 1861">1 : 1</td> <td data-bbox="1094 1776 1214 1861">30"</td> <td data-bbox="1214 1776 1342 1861">60"</td> <td data-bbox="1342 1776 1481 1861">1'30"</td> </tr> <tr> <td data-bbox="352 1861 472 1939">Clinibite Extra Fast</td> <td data-bbox="472 1861 579 1939">95</td> <td data-bbox="579 1861 735 1939">20 µm</td> <td data-bbox="735 1861 890 1939">-0,05%</td> <td data-bbox="890 1861 959 1939">• Light blue • Bleu ciel • Azul claro • Azul claro • Azzurro cielo • Himmelblau • Lichtblauw</td> <td data-bbox="959 1861 1094 1939">1 : 1</td> <td data-bbox="1094 1861 1214 1939">30"</td> <td data-bbox="1214 1861 1342 1939">45"</td> <td data-bbox="1342 1861 1481 1939">1'15"</td> </tr> </tbody> </table>								Technical data Données techniques Datos técnicos Dados técnicos Dati tecnici Technische Daten Technische gegevens	Shore A hardness Dureté Shore A Dureza Shore A Durezza Shore A Shore A Härte Hardheid Shore A	Details reproduction Reproduction des détails Reproducción de los detalles Reprodução de pormenor Riproduzione dei dettagli Detailwiedergabe Reproductie van details	Linear dimensional change Variation dimensionnelle Cambio dimensional Variação dimensional Variazione dimensionale Dimensionsänderung Dimensionale variatie (>24h)	Colour Couleur Color Cor Color Farbe Kleur	Mixing ratio Rapport de mélange Relación de mezcla Relação de mistura Rapporto di miscelazione Anmischverhältnis Mengrapport	Working time * Temps de travail * Tiempo de trabajo * Tempo de trabalho* Tempo de lavorazione* Verarbeitungszeit* Werktijd*	Setting time in mouth Temps de prise en bouge Tiempo de toma en boca Tempo de presa na boca Tempo di permanenza in bocca Verweildauer im mund Tijd mondopname	Total setting time* Temps de prise totale* Tiempo de fraguado total* Tempo de presa total* Tempo di presa totale* Gesamtverarbeitungszeit* Totale opname tijd*	Clinibite Fast	95	20 µm	-0,05%	• Green • Vert • Verde • Grün • Groen	1 : 1	30"	60"	1'30"	Clinibite Extra Fast	95	20 µm	-0,05%	• Light blue • Bleu ciel • Azul claro • Azul claro • Azzurro cielo • Himmelblau • Lichtblauw	1 : 1	30"	45"	1'15"
Technical data Données techniques Datos técnicos Dados técnicos Dati tecnici Technische Daten Technische gegevens	Shore A hardness Dureté Shore A Dureza Shore A Durezza Shore A Shore A Härte Hardheid Shore A	Details reproduction Reproduction des détails Reproducción de los detalles Reprodução de pormenor Riproduzione dei dettagli Detailwiedergabe Reproductie van details	Linear dimensional change Variation dimensionnelle Cambio dimensional Variação dimensional Variazione dimensionale Dimensionsänderung Dimensionale variatie (>24h)	Colour Couleur Color Cor Color Farbe Kleur	Mixing ratio Rapport de mélange Relación de mezcla Relação de mistura Rapporto di miscelazione Anmischverhältnis Mengrapport	Working time * Temps de travail * Tiempo de trabajo * Tempo de trabalho* Tempo de lavorazione* Verarbeitungszeit* Werktijd*	Setting time in mouth Temps de prise en bouge Tiempo de toma en boca Tempo de presa na boca Tempo di permanenza in bocca Verweildauer im mund Tijd mondopname	Total setting time* Temps de prise totale* Tiempo de fraguado total* Tempo de presa total* Tempo di presa totale* Gesamtverarbeitungszeit* Totale opname tijd*																											
Clinibite Fast	95	20 µm	-0,05%	• Green • Vert • Verde • Grün • Groen	1 : 1	30"	60"	1'30"																											
Clinibite Extra Fast	95	20 µm	-0,05%	• Light blue • Bleu ciel • Azul claro • Azul claro • Azzurro cielo • Himmelblau • Lichtblauw	1 : 1	30"	45"	1'15"																											
Storage	<ul style="list-style-type: none"> • 2x50ml + 12 mixing tips. • Storage: temperature between 5° and 27°C. 																																		

Guarantee	3 years
Indication of use	
Legal Information	Class I medical device under European directive 93/42/EEC, including the directives 2007/47/CE CE marking

Version 1

Date: 2020/04/09 par Laetitia DECAUDIN

CFPM - 2 bis chemin du Loup 93290 Tremblay-en-France - FRANCE - Tel. +33 (0)1 45 91 31 82