

ITEM CODE 212.00 Article Mouth mirror handle

Steel	Aspect	Origin	Handle steel	Tip steel	Weight (Kg)
Stainless steel	Shiny	France	AISI303	0	0,027
					<b>Packaging (mm)</b>
					210x30x10

## Marking

Identification		Marking:	
Manufacturer	<i>Prodont-Holliger</i>	Marking:	Laser
Item code		CE	
Lot number			

Cleaning
All means <i>except metallic or abrasive</i> <i>(Ex. Scotch britt)</i>

Sterilisation
all (autoclave, ...) <i>except corrosive products</i> follow cleaning products manufacturers instructions

Warranty
2 years under normal use, except normal wear <i>exchange if defect</i>

Storage
Dry, clean place

Shipping
48 hours on stock

Maintenance
N/A

Medical product surveillance resource
Céline GANDON celine.gandon@acteongroup.com

Last revised version: 2018/01