

Product Specification Sheet - PPR3S



Product Data:

Part Code:	PPR3S
Description:	3 PRICHARD PERIOSTEAL
Instrument Category:	SURGICAL / PERIOSTEALS
Overall Length:	7.00 in (178 mm)
Handle Type:	#522 HEXAGON HANDLE
Weight (Unpkgd/Pkgd):	1,100 oz (31,180 g) / 1,190 oz (33,740 g)
Material Composition:	SEE COMPONENT MATERIALS
	Handle: STAINLESS STEEL Handle Body: N/A Handle Cones: N/A Handle Grips: N/A Working Ends: STAINLESS STEEL
GMDN Term:	Periosteal elevator
GTIN:	10889950020125

Packaging Data:

Material:	Plastic
Style:	Tray with Lid
Weight:	0,090 oz (2,550 g)
Length:	8,020 in (203,708 mm)
Width:	0,730 in (18,542 mm)
Height:	0,400 in (10,033 mm)