

# Fortune Clipper ENTERPRISES

Liaqat Ali Shaheed Road, Street No. 4, Pakka Garha, Sialkot-Pakistan.

## TECHNICAL DATA SHEET

(1)

<b>A. Design Specifications</b>	
1. Size:	As per attached list on Page # 2
2. Name/description of device:	As per attached list on Page # 2
3. Finish:	Mirror/ Satin
4. Markings:	“Stainless” “ Pakistan” “CE” & as per customer demand, if any
5. Intended Use:	Non Sterile Single Use devices to be used in medical
6. Master Samples/illustrations: <b>Note:</b> If customer or any regulatory authority requires a copy of this document, master samples may not be submitted with the document rather illustrations of devices can be submitted.	Available with QA Manager
<b>B. Composition Specifications:</b>	
1. Steel Type used: _	AISI-410
2. Steel Specifications:	<u>ASTM Standard F-899-20</u> Carbon: 0.09 to 0.15 % Manganese 1.0% maximum Phosphorus: 0.040% maximum Sulfur: 0.030% maximum Silicon: 1.0% maximum Chromium: 11.5 to 13.5% Nickle 1.0% maximum
<b>C. Formulation Specifications</b>	
1. Hardness Required	38 TO 43 HRC
<b>D. Component Specifications</b>	
1. Components required	One piece

# Fortune Clipper ENTERPRISES

Liaqat Ali Shaheed Road, Street No. 4, Pakka Garha, Sialkot-Pakistan.

## TECHNICAL DATA SHEET

(2)

### List of Instruments

<b>S. No.</b>	<b>Item Code</b>	<b>Description</b>
1	DEE 20	DENTAL SCALER
	DEE 19	DENTAL SCALER
	CRYER G	DENTAL SCALER
	SINUS 3	DENTAL SCALER
	U 15/30	DENTAL SCALER
	CASTRO DT 18	DENTAL HOLDER
	SINUS 3	DENTAL SCALER
	GRA 3/4	DENTAL SCALER
	GRA 5/6	DENTAL SCALER
	GRA 7/8	DENTAL SCALER
	GRA 9/10	DENTAL SCALER
	GRA 11/12	DENTAL SCALER
	PKT 2	DENTAL SCALER
	EXD 5	DENTAL SCALER
	EXD 7	DENTAL SCALER
	HY 2	DENTAL SCALER
	HY 3	DENTAL SCALER
	LSP 6	DENTAL SCALER
	D 7N	DENTAL SCALER
	D 8N	DENTAL SCALER
	D 9N	DENTAL SCALER
	FOC 0	DENTAL SCALER
	FOC 2	DENTAL SCALER
	LUCAS 2	DENTAL SCALER
	ENDO 3	DENTAL SCALER SILICON
	ENDO 4	DENTAL SCALER SILICON
	ENDO 5	DENTAL SCALER SILICON
	ENDO 6	DENTAL SCALER SILICON
	LECRON	DENTAL SCALER
	GRA 13/14	DENTAL SCALER
	BLD 2	DENTAL SCALER
	HEI TITANE 1	DENTAL SCALER
	HEI TITANE 2	DENTAL SCALER

# Fortune Clipper ENTERPRISES

Liaqat Ali Shaheed Road, Street No. 4, Pakka Garha, Sialkot-Pakistan.

## TECHNICAL DATA SHEET

(3)

### List of Instruments

S. No.	Item Code	Description
	PRICHARD	DENTAL SCALER
	BEALE	DENTAL SCALER
	DEE 20	DENTAL SCALER
	DEE 19	DENTAL SCALER
	CRYER G	DENTAL SCALER
	SINUS 3	DENTAL SCALER
	GRA 3/4	DENTAL SCALER
	GRA 5/6	DENTAL SCALER
	GRA 7/8	DENTAL SCALER
	GRA 9/10	DENTAL SCALER
	GRA 11/12	DENTAL SCALER
	PKT 2	DENTAL SCALER
	EXD 5	DENTAL SCALER
	EXD 7	DENTAL SCALER
	HY 2	DENTAL SCALER
	HY 3	DENTAL SCALER
	LSP 6	DENTAL SCALER
	D 7N	DENTAL SCALER
	D 8N	DENTAL SCALER
	D 9N	DENTAL SCALER
	107	DENTAL PLIER
	FOC 0	DENTAL SCALER
	FOC 2	DENTAL SCALER
	LUCAS 2	DENTAL SCALER
	LECRON	DENTAL SCALER
	GRA 13/14	DENTAL SCALER
	BLD 2	DENTAL SCALER
	HEI TITANE 1	DENTAL SCALER
	HEI TITANE 2	DENTAL SCALER
	PRICHARD	DENTAL SCALER
	BEALE	DENTAL SCALER

# Fortune Clipper ENTERPRISES

Liaqat Ali Shaheed Road, Street No. 4, Pakka Garha, Sialkot-Pakistan.

## TECHNICAL DATA SHEET

(4)

### Material Test Report



**S**tar Steel Casting

G. T. Road, Gujranwala / Pakistan

Tel: 0431-43951 Fax: 0431-45126

Chemical Compositions (%)

Material Supplied by: M/s. Fortune Clipper Enterprises, Pacca Garah, Sialkot.

Dated: 15-10-2020

Batch No: SSC/20/FCE

Size and Article	Weight	C	Si	Mn	P	S	Ni	Cr	Other , If any
Stainless Steel Bar (Grade-410)	1piece	0.10 %	0.41 %	0.42 %	0.023 %	0.024 %	0.45 %	12.66 %	NA

We hereby certify that the above results are true and correct in every detail.

Laboratory Incharge.