


Product	PHOTO RETRACTORS
Reference	Adult : 3011799 Infant : 3011798
Picture	
Description	With an extra-long handle, it is ideal for clinical photography
Sterilization advice	Rinse disinfectants thoroughly before autoclaving. - Autoclavable at 134°C for 18 min. - Cycle : up to 50 x - Do not put metal in contact in the autoclave
Characteristics & Technical data	- Non-reflective - Made of high quality polycarbonate - Sterilizable in autoclave at 134 ° C for 18 min - Available in two sizes: adult and child Store the product preferably in a dry place.
Storage	- Box of 2 adult cheeks retractors - Box of 2 child cheeks retractors
Garantee	More or less 5 years depending on the number of sterilization cycles.
Shell life	n.a
Legal Information	Class I medical device under European directive 93/42/EEC including 2007/47/CE directive CE marking

Version 1.1

Updated 18/09/2018 by LD