## medibase

## TECHNICAL SHEET

Designation	POWDER FREE LATEX EXAMINATION GLOVES SOFT						
Reference	XS -5/6	S-6/7		M - 7/8	Ι	L-8/9	XL-9/10
	3135786	3135787	7 313	5788	31	35789	3135790
Visual product	Inter Examination Gloves Providen Free						
Product description	<ul> <li>Powder-free latex examination gloves, chlorinated.</li> <li>Non-sterile, with micro-textured surface and rolled cuffs.</li> <li>Polymer coated</li> </ul>						
Characteristics and composition	$\begin{tabular}{ c c c c c c c c c c c c c c c c c c c$						110-119
Packaging	Box included 100 gloves (± 2 gloves). 10 boxes per carton.						
Shelf life	5 years from the date of manufacture.						
Indication for use	ISO 374-5-2016 $\overrightarrow{\text{Type B}}$ $\overrightarrow{\text{UPC}}$						
Legal information	Medical Device of class I under MDR 2017/745 Class III personal protective equipment under Regulation 2016/425 CE Marking						

Version 1.0 Update 2021/02/15 by LR