



GC New Metal Strips

Metal Finishing Strips

GC Epitex

Synthetic/Plastic Finishing and Polishing strips



New Metal Strips from GC

Color-coded metal strips, electro-plated with aluminium-oxide particles for approximal contouring, finishing and polishing of all restorative filling materials.

Epitex from GC

Strips for contouring, finishing and polishing in interproximal spaces, suitable for glass ionomers and composites.

GC

GC New Metal Strips

- Advantages**
- › easy to handle, flexible and long-lasting
 - › ultra-thin, yet very strong
 - › can be sterilized or autoclaved repeatedly
 - › various grades and widths (2,6-3,3-4,0 mm)

- Indications**
- › N° 200 (red) - preparation
 - › N° 300 (blue) - final preparation
 - › N° 600 (green) - finishing
 - › N° 1.000 (yellow) - polishing



GC EpiteX

- Advantages**
- › Extremely thin and flexible material
 - › Excellent adaptability even between tight interproximal contacts
 - › Very tear resistant
 - › Abrasive particles are firmly embedded into the strip surface
 - › Four grades
 - › Handy strip dispenser
 - › Also available as matrix strips

- Indications**
- › Coarse contouring (blue)
 - › Medium finishing (green)
 - › Fine polishing (grey)
 - › Extra fine polishing (pink)



Assortment pack

GC New Metal Strips	
Assortment pack	
000267	GC New Metal Strips assortment Content 12 pcs. (1 pcs. of each width and grain)
Content: 12 pcs. (4pcs. each of widths 2.6mm, 3.3mm, 4.0mm)	
000251	GC New Metal Strips Red Mesh N° 200 (140µm Rough Preparation)
000255	000255 GC New Metal Strips Blue Mesh N° 300 (90µm Preparation)
000259	000259 GC New Metal Strips Green Mesh N° 600 (70µm Finishing)
000263	000263 GC New Metal Strips Yellow Mesh N° 1000 (50µm Polishing)
6mm Refill packs: (12pcs. of 6mm wide)	
000268	GC New Metal Strips Red Mesh N° 200 (140µm Rough Preparation)
000269	GC New Metal Strips Blue Mesh N° 300 (90µm Preparation)
000270	GC New Metal Strips Green Mesh N° 600 (70µm Finishing)
000271	GC New Metal Strips Yellow Mesh N° 1000 (50µm Polishing)
GC EpiteX	
Assortment pack	
000404	GC EpiteX, Starter Kit
Content: Stand with dispenser containing all available grains, 10m each (coarse-blue, medium-green, fine-grey, extra fine-red)	
Unit refill packs (Content: 10m)	
000409	GC EpiteX Translucent matrix, 10 m
000405	GC EpiteX Coarse (blue), 10 m
000407	GC EpiteX Medium (green) 10 m
000406	GC EpiteX Fine (grey), 10 m
000408	GC EpiteX Extra fine (red), 10 m

GC EUROPE N.V.

Head Office
 Researchpark Haasrode-Leuven 1240
 Interleuvenlaan 33, B-3001 Leuven
 Tel. +32.16.74.10.00, Fax. +32.16.40.48.32
 info.gce@gc.dental
<https://europe.gc.dental>

GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell
 Buckinghamshire MK16 8JS
 United Kingdom
 Tel. +44.1908.218.999, Fax. +44.1908.218.900
 info.uk@gc.dental
<https://europe.gc.dental/en-GB>

