




Technical sheet

| Designation | INTERDENTAL WOODEN WEDGES | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|----------------------------|--|------------------------|------------------|---------|---------|--------------|-------------|------------------------|------------------|---------------------|----|-----|-----|----------------|----|-----|-----|---------------|----|-----|-----|---------------|----|-----|-----|-----------------|----|-----|-----|
| Reference | Turquoise | Green | Orange | Yellow | White | | | | | | | | | | | | | | | | | | | | | | | | |
| | XS | S | M | L | XL | | | | | | | | | | | | | | | | | | | | | | | | |
| | 0286866 | 0286867 | 0124214 | 0124215 | 3021163 | | | | | | | | | | | | | | | | | | | | | | | | |
| Product description | Boxwood interdental wedges for ideal tooth-matrix adaptation. | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Picture |  | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Characteristics | <ul style="list-style-type: none"> - No sliding after insertion. - Anatomical form guarantees successful contact. | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Indications | <table border="1" data-bbox="523 1267 1382 1789"> <thead> <tr> <th data-bbox="523 1267 743 1435">Wooden wedge</th> <th data-bbox="743 1267 1007 1435">Length (mm)</th> <th data-bbox="1007 1267 1177 1435">Base width ϕ (mm)</th> <th data-bbox="1177 1267 1382 1435">Base height (mm)</th> </tr> </thead> <tbody> <tr> <td data-bbox="523 1435 743 1509">X-Small (Turquoise)</td> <td data-bbox="743 1435 1007 1509">11</td> <td data-bbox="1007 1435 1177 1509">1.2</td> <td data-bbox="1177 1435 1382 1509">1.6</td> </tr> <tr> <td data-bbox="523 1509 743 1583">Small (Orange)</td> <td data-bbox="743 1509 1007 1583">12</td> <td data-bbox="1007 1509 1177 1583">1.6</td> <td data-bbox="1177 1509 1382 1583">1.6</td> </tr> <tr> <td data-bbox="523 1583 743 1657">Medium (Vert)</td> <td data-bbox="743 1583 1007 1657">13</td> <td data-bbox="1007 1583 1177 1657">2.0</td> <td data-bbox="1177 1583 1382 1657">2.0</td> </tr> <tr> <td data-bbox="523 1657 743 1731">Large (Jaune)</td> <td data-bbox="743 1657 1007 1731">15</td> <td data-bbox="1007 1657 1177 1731">2.0</td> <td data-bbox="1177 1657 1382 1731">2.2</td> </tr> <tr> <td data-bbox="523 1731 743 1792">X-Large (Blanc)</td> <td data-bbox="743 1731 1007 1792">17</td> <td data-bbox="1007 1731 1177 1792">2.3</td> <td data-bbox="1177 1731 1382 1792">2.3</td> </tr> </tbody> </table> | | | | | Wooden wedge | Length (mm) | Base width ϕ (mm) | Base height (mm) | X-Small (Turquoise) | 11 | 1.2 | 1.6 | Small (Orange) | 12 | 1.6 | 1.6 | Medium (Vert) | 13 | 2.0 | 2.0 | Large (Jaune) | 15 | 2.0 | 2.2 | X-Large (Blanc) | 17 | 2.3 | 2.3 |
| Wooden wedge | Length (mm) | Base width ϕ (mm) | Base height (mm) | | | | | | | | | | | | | | | | | | | | | | | | | | |
| X-Small (Turquoise) | 11 | 1.2 | 1.6 | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Small (Orange) | 12 | 1.6 | 1.6 | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Medium (Vert) | 13 | 2.0 | 2.0 | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Large (Jaune) | 15 | 2.0 | 2.2 | | | | | | | | | | | | | | | | | | | | | | | | | | |
| X-Large (Blanc) | 17 | 2.3 | 2.3 | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Composition | - Wood | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Packaging | - Kit of 100 wooden wedges | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Legal information | Medical Device of class I under European Directive 93/42 /EEC including Directive 2007/47/CE CE Marking | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Version 2.0

Updated 2016/06/28 by Scarlet Ceausoiu

CFPM - 2 bis chemin du Loup 93290 Tremblay-en-France - FRANCE - Tel. +33 (0)1 45 91 31 82