

LARIDENT S.R.L.UNIPERSONALE
VIA LAMANEIGRA, 12 R
161030 TRIBOGNA (GE) ITALIA
TEL.0185938068 FAX.0185938595 info@larident.it

TECNICAL SHEET

OF PRODUCT

PRODUCT NAME	ORTHODONTIC BOX WITH VENTILATION HOLES AND A HOLE FOR THE CORD AVAILABLE IN SEVERAL COLOURS
CODE	G3
MATERIAL	POLYPROPYLENE-HOMOPOLYMER AND COLORED MASTER
THERMAL AND PHYSICAL PROPERTY	AUTOCLAVABLE UP TO 140°C. ANTIBREAKING ZIP MEASUREMENTS: CM 8*7 H4 NO PHTHALATES IT RESISTS TO ALL DISINFECION LIQUIDS
ACCESSORY	G20 COLORED CORD FOR G3
PACKING	IT IS SOLD IN PACKING OF 10 PCS, IN A PLASTIC BAG AVAILABLE IN SEVERAL COLOURS : WHITE – RED –BLUE-GREEN-FUSCIA- YELLOW – YELLOW FLOU- VIOLET- ORANGE – LIGHT BLUE AND SPECIAL GLITTER COLOURS
SAFETY DATA SHEET	N/A
WARNING	NONE
ANALYSIS OF RISK R	N/A
CLASSIF. CEE 93/42	NOT A MEDICAL DEVICE

N/A NO APPLICABLE