

▼ Mandibular		Torq.	Ang.	Ancho	R/L	.022" Slot
$\frac{1}{2} \frac{1}{2}$	Anterior Inferior	-6°	0°	.096"	U	N-753
$\frac{3}{3}$	Cúspide con gancho	-6°	+3°	.119"	R	N-745
					L	N-749
$\frac{3}{3}$	Cúspide con gancho	0°	+3°	.119"	R	N-746
					L	N-750
$\frac{4}{4}$	1 st Bicúspide c/gancho y off	-12°	+2°	.119"	R	N-747
					L	N-751
$\frac{5}{5}$	2 nd Bicúspide c/gancho y off	-17°	+2°	.119"	R	N-748
					L	N-752

▲ Maxilar		Torq.	Ang.	Ancho	R/L	.022" Slot
$\frac{1}{1}$	Central	+17°	+4°	.130"	R	N-733
					L	N-740
$\frac{2}{2}$	Lateral	+10°	+8°	.109"	R	N-734
					L	N-741
$\frac{3}{3}$	Cúspide con gancho	-7°	+8°	.117"	R	N-735
					L	N-742
$\frac{3}{3}$	Cúspide con gancho	0°	+8°	.117"	R	N-736
					L	N-743
$\frac{4}{5} \frac{4}{5}$	Bicúspide c/gancho y off	-7°	0°	.120"	R	N-737
					L	N-744