



Technical sheet

Désignation	TEMPORARY CROWN & BRIDGE										
References		SHADE	REFERENCE								
		A2	3124456								
		A3	3124457								
		A3,5	3124458								
		BLEACH	3143475								
Product description	<ul style="list-style-type: none"> Self-curing resin for the preparation and manufacture of temporary prostheses (crowns, bridges, inlays/onlays, veneers) 										
Picture											
Characteristics and technical data	<ul style="list-style-type: none"> Packaging : a 98g cartridge + 10 mixing tips 3 available colours: A2, A3, A3.5 Compatible mixing gun: type 1:1 										
	<table border="1"> <thead> <tr> <th>Item</th> <th>Water sorption ($\mu\text{m}/\text{mm}^2$)</th> <th>Solubility ($\mu\text{m}/\text{mm}^2$)</th> <th>Flexural strength (MPa)</th> </tr> </thead> <tbody> <tr> <td>Test result</td> <td>>40</td> <td>>75</td> <td><50</td> </tr> </tbody> </table>				Item	Water sorption ($\mu\text{m}/\text{mm}^2$)	Solubility ($\mu\text{m}/\text{mm}^2$)	Flexural strength (MPa)	Test result	>40	>75
Item	Water sorption ($\mu\text{m}/\text{mm}^2$)	Solubility ($\mu\text{m}/\text{mm}^2$)	Flexural strength (MPa)								
Test result	>40	>75	<50								
Advantages	<ul style="list-style-type: none"> Easy and fast handling High compressive and bending strength Easy to finish and polish Low polymerization shrinkage Very aesthetic final result Several shades available Standard mixing gun 										
Legal information	Class IIa medical device.										

Version 4.0
Updated 14/11/2022 by RG

CADENCE - 2 bis chemin du Loup 93290 Tremblay-en-France - FRANCE - Tel. +33 (0)1 45 91 31 82