



Technical Datasheet

Product	SUPRAPOWDER					
References	<table border="1"> <thead> <tr> <th>MENTHE</th> <th>CITRON</th> </tr> </thead> <tbody> <tr> <td>3140107</td> <td>3140106</td> </tr> </tbody> </table>	MENTHE	CITRON	3140107	3140106	
MENTHE	CITRON					
3140107	3140106					
Description	Prophylactic powder based on sodium bicarbonate for the supragingival removal of stubborn deposits and stains on tooth enamel.					
Visual						
Characteristics	<ul style="list-style-type: none"> – Made of sodium bicarbonate (particle size = 40 µm). – Supragingival removal of stubborn deposits and stains on tooth enamel. – Compatible with all air polishers. – Soluble in water. 					
Composition	Sodium hydrogen carbonate based					
Packaging	Bottle of 300g. 2 flavors available: mint and lemon.					
Storage	Storage temperature: minimum -10°C / maximum 30 °C					
Shelf life	42 months					
Legal Information	Medical Device under Class I, 2017/745 – CE marking					