



MICRO-TOUCH[®] Coated

Latex, powder-free, examination glove coated for quick and frequent donning

GENERAL DESCRIPTION	
Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	White
Shape	Ambidextrous
Internal Glove Surface	Polyacrylate Coated
External Glove Surface	Textured Fingers
Cuff Style	Beaded

PHYSICAL PROPERTIES			
	Finger	Palm	Cuff
Thickness Typical Average (mm/mils)	0.14 / 5.5	0.11 / 4.3	0.09 / 3.5
Glove Length Typical (mm/in)	240 / 9.5		
Strength Minimum	Before Aging / Fresh		
	ASTM ¹	EN ²	
Elongation At Break (%)	≥650	-	
Force At Break (N)	-	≥ 6	
Tensile Strength (MPa)	≥18	-	
	After Aging / After Challenge		
	ASTM ¹	EN ²	
Elongation At Break (%)	≥500	-	
Force At Break (N)	-	≥ 6	
Tensile Strength (MPa)	≥14	-	

ALLERGENIC SAFETY PROFILE	
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)
Primary Skin Irritation	Considered a negligible irritant as per ISO 10993-10
Skin Sensitization	No evidence of delayed dermal contact sensitization as per ISO 10993-10
Protein Level	30 µg/g or less of total extractable protein

SAFETY STANDARDS	
Freedom From Holes (AQL)	1.5 AQL
Sterilization Method	Not Applicable
Tested For Use With Chemotherapy Drugs	No
Viral Penetration	Passes ASTM F1671 using Phi X 174, Passes ISO 16604 using Phi X 174

STANDARDS AND CERTIFICATIONS	
Product Standards	CE 0493, Category III, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 11193-1:2008, ISO 13485, ISO 10993-10:2010
Quality / Environmental Standards	ISO 13485, ISO 9001
Product Certification	Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks

Contact your Ansell representative for ordering or more information.

North America

US: 800 952 9916

CA: 1-844-494-7854

insidesalesus@ansell.com

Central & South America

+52(442) 296 20 50

cslac@ansell.com

Brazil

+55-11-3356-3100

luvas.medicas@ansell.com

Europe, Middle East & Africa

+32 (0)2 528 74 00

info@ansell.eu

Asia Pacific

+603 8310 6688

apac.medical@ansell.com

Australia & New Zealand

+61 3 9270 7270

protection@ap.ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. © 2022 Ansell Limited. All Rights Reserved.

MICRO-TOUCH[®] Coated

Latex, powder-free, examination glove coated for quick and frequent donning

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case *80 gloves per dispenser (Size XL) *800 gloves per case (Size XL)
Shelf Life	3 Years
Shared Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

ANSELL PRODUCT CODES

Size	Code	Size	Code
XS (5.5 - 6)	553301	S (6.5 - 7)	553302
M (7.5 - 8)	553303	L (8.5 - 9)	553304
XL (9.5 - 10)	553305		

CE 0493



EN ISO 374-1:2016
Type B



EN ISO 374-5



ISO 9001



ISO 13485

Contact your Ansell representative for ordering or more information.

North America

US: 800 952 9916

CA: 1-844-494-7854

insidesalesus@ansell.com

Central & South America

+52(442) 296 20 50

csllac@ansell.com

Brazil

+55-11-3356-3100

luvas.medicas@ansell.com

Europe, Middle East & Africa

+32 (0)2 528 74 00

info@ansell.eu

Asia Pacific

+603 8310 6688

apac.medical@ansell.com

Australia & New Zealand

+61 3 9270 7270

protection@ap.ansell.com

ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. © 2022 Ansell Limited. All Rights Reserved.